

FALL PREVENTION ALARMS

This alarm is to keep an eye on the patient who is on a wheelchair or a sick bed. Once the patient departs from the chair or bed pad, the alarm will be activated for calling nurse to help. This device can also be connected to nursing station interface.

這是永遠不會疲勞打瞌睡的看護，只要看護的病人從輪椅或病床上離開，立即發出警告聲求援跌落告知警報器，並可以連接護理站呼叫系統介面，轉知護理人員

The plug connects to the nursing station interface



• LD-166 Patent No. M273054



• PD-CHAIR Patent No. 195100



• PD-BED Patent No. 195100


Special anti-slippery bag designed for attaching to any part of the bed frame



• LD-166


Reset button

Inside the lockable sliding cover, there are three switches for adjusting volume, tempo & tone if required. The protective cover will prevent the wires inside from disconnecting



• PD-MAT

Safety floor pad alarm: preventing incoherent patient from leaving their bed. Once the patient steps on the pad, the alarm will continuously alert the nurses for assistance



Care-alert Alarm and Holder

Activated when patient steps on