

## FALL PREVENTION ALARMS

The most simple fall prevention alarm that can be used for retirement communities as well as family homes.

最簡潔有效看護警報裝置，適用於自行看護場合或家庭內



Use a locking-clip on the end of pulling strap to fix the alarm firmly on clothes.

Separating to activate



Once the patient falls down from the bed, separating the magnet switch, the alarm will be activated.



### • LD-68PS or LD-166

The reset button will act as a deactivator after the alarm being triggered. It replaces the need of on and off switch and eliminates the risk of forgetting to turn it on.

## BACKPACK ALARM

Once the alarm is activated, a very loud siren is bursting out from both sides of the backpack and very bright LED strobe is flashing at the middle of the lower portion of the backpack. The resonance of the siren from these twin piezos will be combined together to make it louder after few seconds of continuous activation.

背包警報器：遇狀況時，隨意拉脫啟動插梢，震耳欲聲的警報聲由背包兩旁發出又由於共震的加強效果，警報音會越來越震撼同時強烈紅色警示燈閃爍指示所在位置（尤其於夜間），適用於登山，夜歸婦女



### • BP-1 Patent No. 179783

### Key Specifications

- Tone : Combination of twin piezo sirens' resonance
- Siren output : 130 dB ( at 1 M )
- Stand-by current consumption : 1.0uA
- Battery : " AA " X 8